

Dandle WRAP™ and Dandle WRAP Stretch



INTENDED USE / CLINICAL BENEFIT:

The Dandle WRAP and Dandle WRAP *Stretch* are reusable positioning aids that can be used with hospitalized infants to provide neurodevelopmental support while the baby is sleeping, eating or being bathed or held. It can also be used to provide support during procedures or while the baby is being handled. The Dandle WRAP provides additional support for low tone and premature infants while adhering to the Back to Sleep guidelines.

WARNINGS / CONTRAINDICATIONS:

- **Hospital use only**
- **Do not place wings or pouch around baby's neck or face or use in any manner that might restrict infant's breathing**
- **Between patients, wash according to hospital policy**
- **Monitor infant's temperature to prevent overheating**
- **Be sure Velcro® is attached securely by pressing it together from both sides. Avoid contact between Velcro hook and baby's skin**

CLEANING INSTRUCTIONS:

- Machine wash in water temperatures up to 71°C. Use only non-chlorine bleach. The product will hold up better if all Velcro fasteners are attached during laundering. Tumble dry using medium heat.



INSTRUCTIONS FOR USE:

Select the correct size Dandle WRAP according to the infant's size and body proportions.

1. Place Dandle WRAP on a supportive surface with the wings out to the sides and the pouch opened up at the bottom. The Velcro loop should be on the outside of the back.
2. Place baby supine on the WRAP with the top of the shoulders lined up just below the top of the WRAP.
3. Bring upper extremities into midline flexion with hands near face. Pull the shorter wrap down over the shoulder and across the arms, rounding shoulder forward. Pull the longer wrap over top, rounding shoulders and containing arms comfortably. Fasten to the Velcro square on the backside of the positioner. Adjust according to baby's needs and medical condition.
4. Flex baby's lower extremities into a neutral midline position, supporting the pelvis in a posterior tilt.
5. Safely arrange any tubing and medical equipment so they may exit the side.
6. Bring pouch up over the legs and feet. Fasten to the backside. Adjust according to baby's needs and medical condition.
7. May be used in any position (prone, supine or side-lying) with medical supervision.

STORAGE AND DISPOSAL:

No special requirements.

BACK TO SLEEP GUIDELINES for use with Dandle WRAP and Dandle WRAP *Stretch*:

Follow all instructions except:

- a. Step 3: Ensure that WRAP does not become loose around baby's face by either securing it snugly with the fasteners or placing the shorter wrap underneath the baby's armpits rather than over the shoulders to ensure baby will not nestle down into the WRAP.
- b. Step 7: Baby should always be placed on his/her back to sleep.

DEHP-free, BPA-free, not made with any natural rubber latex.

Please contact Dandle•LION Medical with any questions or comments or to report any issue with the use of our product.

Rev6:0821