



# Dandle WRAP™ Stretch

## Dynamic Positioning Aid

The Dandle WRAP *Stretch* offers many swaddling benefits over traditional blankets:

- High performance fabric with 3 way stretch allows freedom of movement with return to midline flexion
- A pouch to maintain lower extremity flexion
- Wide wings that support the shoulders in flexion and encourage hand to mouth self-comforting behaviors
- Velcro® fasteners on pouch and wings that can be adjusted independently of each other to meet baby's individual needs
- Fabric that breathes, yet is absorbent, making it ideal for agitated babies
- Plenty of room between the pouch and the wings for lines or tubing
- Meets AAP Back to Sleep guidelines

Babies love to be swaddled! By containing motor activity and mimicking the security of a mother's womb, swaddling comforts babies and results in less crying, more stable heart rates and the ability to sleep through the startle reflex. In a systematic review of swaddling published in *Pediatrics* in 2007, authors found "swaddled infants arouse less and sleep longer. Preterm infants have shown improved neuromuscular development, less physiologic distress, better motor organization, and more self-regulatory ability when they are swaddled... (Swaddling) is supportive in cases of neonatal abstinence syndrome and infants with neonatal cerebral lesions."

Swaddling: A Systematic Review van Sleuwen et al, *Pediatrics*, Vol. 120, No. 4 October 1, 2007, pp. e1097-1106.

  
Dandle•LION™  
MEDICAL



# Dandle WRAP Stretch



Dandle WRAP Stretch is made with fabulous high performance fabric!


- 3 way stretch that returns to its original shape
- Soft, smooth, breathable, wickable, and highly absorbent
- Does not insulate body heat, making it ideal for use in hot and humid environments
- Washes well even in industrial machines with minimal shrinkage and pilling



- **Helps During Feeding** – keeps baby organized and contained. Staff can feel respirations through the lightweight fabric and easily see the head and neck in order to visualize sucking and swallowing



- **Calms Even the Most Irritable Infants** – adjustable swaddling tension and secure Velcro help to calm babies and maintain a flexed and midline position



- **Supports During Procedures and Diaper Changes** – wings and pouch can be adjusted independently of each other to keep baby calm and support baby's needs

D 14531 B12	Dandle WRAP Stretch – XS: <1000 g.	12/case
D 14532 B12	Dandle WRAP Stretch – S: 1000-1800 g.	12/case
D 14533 B12	Dandle WRAP Stretch – M: 1800-2500 g.	12/case
D 14534 B12	Dandle WRAP Stretch – L: 2500-3600 g.	12/case
D 14535 B12	Dandle WRAP Stretch – XL: 3600-4800 g.	12/case
D 14536 B12	Dandle WRAP Stretch – XXL: 4800 g.+	12/case