

Available in
two sizes!



48" Long

40" Long



Dandle•LION Cozy Cub™

Bead Boundary in 2 Sizes

The Cozy Cub:

- Provides containment while allowing full visualization, making it ideal for the ECMO, cardiac or surgical infants
- Allows immediate access to the infant in emergencies
- Useful for patients with unique or critical positioning needs
- Easily adaptable for use with medical equipment
- Capable of providing calming flexion during transport
- Sizes:
 - L: 48" long, 6 lbs
 - S: 40" long, 3.5 lbs
- Packaged with a disposable cover
- Additional printed cloth covers and disposable white covers available separately.

The Cozy Cub is a pellet-filled boundary that comes in 2 sizes to help meet the needs of all hospitalized babies. The small internal pellets are smooth and almost feel silky. Unlike more formed boundaries, the Cozy Cub can get into the nooks and crannies to provide support exactly where it is needed. Developed originally in the larger size for ECMO and post-surgical patients, both the large and small Cozy Cub can be used for preemies, agitated babies of all sizes, or any baby who just needs a little 'more.'



Loop at one end helps caregiver adjust device to fit any size baby

Ordering Information

D 17205 B4 | Cozy Cub w/white non-woven cover - Small | Qty: 4/case

D 17005 B4 | Cozy Cub w/white non-woven cover - Large | Qty: 4/case

D 17215C B25 | Printed cotton covers - Small | Qty: 25/box

D 17015C B25 | Printed cotton covers - Large | Qty: 25/box

D 17225C B25 | Disposable white non-woven covers - Small | Qty: 25/box

D 17025C B25 | Disposable white non-woven covers - Large | Qty: 25/box

Made in the USA

Dandle•LION Medical | 22 Shelter Rock Lane | Danbury, CT 06810 | 203-791-9001 | www.dandleLIONmedical.com

Rev8:0322

Dandle•LION™
MEDICAL