

Polly PAL™

Infant Positioning Aid

Not made with natural rubber latex. DEHP free.


Polly PAL™ is a flexible positioning aid that can be used in a variety of ways. It is filled with plastic micro discs, that feel silky to the touch.

WARNINGS:

- Hospital use only.
- Special caution should be taken to avoid placing the full weight or too many micro discs on top of small infants.
- Do not place near the baby's neck or face or use in any manner that might restrict the baby's airway.
- To prevent cross-contamination between patients, wash and dry according to hospital policies.
- If micro discs should ever become loose or if a hole occurs in the fabric cover, discard immediately.

CLEANING INSTRUCTIONS:

The Polly PAL is washable and reusable product.



INSTRUCTIONS FOR USE:

The Polly PAL is a flexible positioning aid that can be used in a variety of ways.

- Place around the baby's head or feet to create a boundary, help to maintain a position or give a feeling of security.
- Use it to prop up a soother or support tubing.
- Move the micro discs into different places to make parts or sections fuller or less full to accommodate equipment or the baby's positioning needs.

DISPOSAL and STORAGE:

No special requirements.

Please contact DandleLION medical with any questions or comments or to report any issue with the use of our product.

Rev1:0621

Polly PAL™

Infant Positioning Aid

Not made with natural rubber latex. DEHP free.


Polly PAL™ is a flexible positioning aid that can be used in a variety of ways. It is filled with plastic micro discs, that feel silky to the touch.

WARNINGS:

- Hospital use only.
- Special caution should be taken to avoid placing the full weight or too many micro discs on top of small infants.
- Do not place near the baby's neck or face or use in any manner that might restrict the baby's airway.
- To prevent cross-contamination between patients, wash and dry according to hospital policies.
- If micro discs should ever become loose or if a hole occurs in the fabric cover, discard immediately.

CLEANING INSTRUCTIONS:

The Polly PAL is washable and reusable product.



INSTRUCTIONS FOR USE:

The Polly PAL is a flexible positioning aid that can be used in a variety of ways.

- Place around the baby's head or feet to create a boundary, help to maintain a position or give a feeling of security.
- Use it to prop up a soother or support tubing.
- Move the micro discs into different places to make parts or sections fuller or less full to accommodate equipment or the baby's positioning needs.

DISPOSAL and STORAGE:

No special requirements.

Please contact DandleLION medical with any questions or comments or to report any issue with the use of our product.

Rev1:0621