



Product Catalog



Table of Contents

Developmental Positioning Aids

Dandle® ROO2	3
Dandle® ROO Lite	4
DandleLION® Snuggly ROO	5
Dandle® WRAP Stretch	6
Crown Cradle™	7
Cozy Cub™	8
Dandle® PAL, Polly PAL, & Mini PAL™	9
DandleLION® Gel Pillows	10
Dandy Prone Pad™	11

Organic Skin Care

Bubbly Wash	12
Diaper Balm	12
Nurturing Oil	12
Silky Cream	12

Swaddled Bathing

DandleLION® Bathing System	13
----------------------------	----

Baby Care

Tranquilo® Soothing Mat	14
Thermal Hats	14
DandleLION® Sweet Cheeks	15
DandleLION® KISSES	15

Contact Us

Contact Information	20
---------------------	----

Dandle® ROO2

A developmental positioning aid



A developmental positioning aid made for hospitalized infants made with breathable, wickable, 3-way stretch material. Comes with base unit (headroll, wings and pouch) and 1 multi-purpose roll. Accessories to meet each baby's unique positioning and developmental needs.

SKU	DESCRIPTION	QTY
D 14510 B12	Dandle Roo2 - XS: <1000 g.	12/case
D 14510 ROLL B25	Dandle Roo2 - XS: <1000 g. Multi-purpose rolls	25/case
D 14510 BRM B25	Dandle Roo2 - XS: <1000 g. Detachable brim	25/case
D 14510 GP B12	Dandle Roo2 - XS: Gel Pillow w/ printed cotton cover - fits 6" x 9" pillows	12/case
D 14510 GPC B25	Dandle Roo2 - XS: Gel Pillow printed cotton covers - fits 6" x 9" pillows	25/case
D 14520 B12	Dandle Roo2 - S: 1000-1800 g.	12/case
D 14520 ROLL B25	Dandle Roo2 - S: 1000-1800 g. Multi-purpose rolls	25/case
D 14520 BRM B25	Dandle Roo2 - S: 1000-1800 g. Detachable brim	25/case
D 14520 GP B12	Dandle Roo2 - S: Gel Pillow w/ printed cotton cover - fits 6" x 11" pillows	12/case
D 14520 GPC B25	Dandle Roo2 - S: Gel Pillow printed cotton covers - fits 6" x 11" pillows	25/case
D 14530 B12	Dandle Roo2 - M: 1800-2500 g.	12/case
D 14530 ROLL B25	Dandle Roo2 - M: 1800-2500 g. Multi-purpose rolls	25/case
D 14530 BRM B25	Dandle Roo2 - M: 1800-2500 g. Detachable brim	25/case
D 14530 GP B12	Dandle Roo2 - M: Gel Pillow w/ printed cotton cover - fits 6" x 12" pillows	12/case
D 14530 GPC B25	Dandle Roo2 - M: Gel Pillow printed cotton covers - fits 6" x 12" pillows	25/case
D 14550 B12	Dandle Roo2 - Starter Pack - 4 each size: XS, S, M	12/case
D 0981 B25	DandleLION laundry bags - 12" x 15" blue mesh	25/case

Dandle® ROO Lite

A disposable developmental positioning aid



The Dandle ROO Lite may be used during phototherapy or when a disposable product better meets the needs of the baby. The ROO Lite may be hand laundered. 73% of phototherapy light passes through 1 layer of the 2 oz. cotton wings and pouch.

SKU	DESCRIPTION	QTY
D 11031 B12	Dandle ROO Lite - XS: <1000 g.	12/case
D 11031 BRM B25	Dandle ROO Lite Brim - XS: <1000 g.	25/case
D 0951 ROLL B25	Dandle ROO - XS: multi-purpose rolls (for ROO or ROO Lite)	25/case
D 0951 GP B12	Dandle ROO - XS: Gel Pillow w/ White cotton cover (for ROO or ROO Lite) - 6" x 9"	12/case
D 0951 GPC B25	Dandle ROO - XS: Gel Pillow covers - White cotton	25/case
D 12551 GPC B25	Dandle ROO - XS: Disposable Gel Pillow covers	25/case
D 11032 B12	Dandle ROO Lite - S: 1000-1800 g.	12/case
D 11032 BRM B25	Dandle ROO Lite Brim - S: 1000-1800 g.	25/case
D 0952 ROLL B25	Dandle ROO - S: multi-purpose rolls (for ROO or ROO Lite)	25/case
D 0952 GP B12	Dandle ROO - S: Gel Pillow w/ White cotton cover (for ROO or ROO Lite) - 6" x 11"	12/case
D 0952 GPC B25	Dandle ROO - S: Gel Pillow covers - White cotton	25/case
D 12552 GPC B25	Dandle ROO - S: Disposable Gel Pillow covers	25/case
D 11033 B12	Dandle ROO Lite - M: 1800-2500 g.	12/case
D 11033 BRM B25	Dandle ROO Lite Brim - M: 1800-2500 g.	25/case
D 0953 ROLL B25	Dandle ROO - M: multi-purpose rolls (for ROO or ROO Lite)	25/case
D 0953 GP B12	Dandle ROO - M: Gel Pillow w/ White cotton cover (for ROO or ROO Lite) - 6" x 12"	12/case
D 0953 GPC B25	Dandle ROO - M: Gel Pillow covers - White cotton	25/case
D 12553 GPC B25	Dandle ROO - M: Disposable Gel Pillow covers	25/case
D 11034 B12	Dandle ROO Lite - L: 2500-3600 g.	12/case
D 11034 BRM B25	Dandle ROO Lite Brim - L: 2500-3600 g.	25/case
D 0981 B25	DandleLION laundry bags - 12" x 15" blue mesh	25/case

DandleLION® Snuggly ROO

Cloth and Disposable Positioning Aids



The Snuggly ROO is a cozy positioning nest with a padded footroll and two straps that keep babies flexed, contained and comfortable with their extremities in midline. The disposable Snuggly ROO is ideal for babies at admission to keep them supported during procedures, and also a great alternative to traditional cloth positioning aids when laundry services are not available. The cloth Snuggly ROO can be washed and reused between patients and is useful during times when greater visualization of the infant is important.

SKU	DESCRIPTION	QTY
D 21621 B12	DandleLION Snuggly ROO - XS:<1000 g.	12/case
D 21622 B12	DandleLION Snuggly ROO - S: 1000-1800 g.	12/case
D 21623 B12	DandleLION Snuggly ROO - M: 1800-2500 g.	12/case
D 21624 B12	DandleLION Snuggly ROO - L: 2500-3600 g.	12/case
D 21031 B12	DandleLION Snuggly ROO Disposable - XS:<1000 g.	12/case
D 21032 B12	DandleLION Snuggly ROO Disposable - S: 1000-1800 g.	12/case
D 21033 B12	DandleLION Snuggly ROO Disposable- M: 1800-2500 g.	12/case
D 21034 B12	DandleLION Snuggly ROO Disposable - L: 2500-3600 g.	12/case

Dandle® WRAP Stretch

A dynamic positioning aid



A one-piece positioning aid that supports developmental positioning and provides overall containment for the hospitalized baby. Breathable, wickable, 3-way stretch material holds baby in flexion and provides excellent containment for agitated babies.

SKU	DESCRIPTION	QTY
D 14531 B12	Dandle WRAP Stretch - XS: >1000 g.	12/case
D 14532 B12	Dandle WRAP Stretch - S: 1000-1800 g.	12/case
D 14533 B12	Dandle WRAP Stretch - M: 1800-2500 g.	12/case
D 14534 B12	Dandle WRAP Stretch - L: 2500-3600 g.	12/case
D 14535 B12	Dandle WRAP Stretch - XL: >3600-4800 g.	12/case
D 14536 B12	Dandle WRAP Stretch - XXL: 4800+ g.	12/case
D 14538 B12	Dandle WRAP Stretch Sample Pack: 4 ea. M, L, XL	12/case
D 0981 B25	DandleLION laundry bags - 12" x 15" blue mesh	25/case

Crown Cradle™

Preemie orthotic device to support development of infant's natural head shape



The Crown Cradle may reduce the potential impact to the infant's head by redistributing contact pressure, supporting the cranium, and providing a pathway for more symmetrical skull growth. It helps the bedside caregiver to confidently developmentally position the baby, provide appropriate alignment of the spine, and support the developing brain and skull.

SKU	DESCRIPTION	QTY
D 18510 B2	Crown Cradle System SMALL (Infants <1800 g), Includes Mattress, Mattress Plug, Mattress Cover, Foam Insert Sizes 1-4	2/case
D18512 B2	Crown Cradle Replacement Mattresses, SMALL, Includes Mattress Plugs	2/case
D18513 B6	Crown Cradle Replacement Mattress PLUG, SMALL	6/case
D18516 B6	Crown Cradle Replacement Mattress COVERS, SMALL	6/case
D 18525 PK	Crown Cradle Replacement Foam Inserts, SMALL Mattress Multi Size Pack, Size 1 - 4, ONE of each size	4/pack; 1 pack/case
D 18530 B2	Crown Cradle System LARGE (Infants >1800 g), Includes Mattress, Mattress Plug, Mattress Cover, Foam Insert Sizes 5-8	2/case
D 18532 B2	Crown Cradle Replacement Mattresses, LARGE, Includes Mattress Plugs	2/case
D 18533 B6	Crown Cradle Replacement Mattress PLUG, LARGE	6/case
D 18536 B6	Crown Cradle Replacement Mattress COVERS, LARGE	6/case
D 18545 PK	Crown Cradle Replacement Foam Inserts, LARGE Mattress Multi Size Pack, Size 5 - 8, ONE of each size.	4/pack; 1 pack/case

Cozy Cub™

A bead boundary in 3 sizes



The Cozy Cub is a pellet-filled boundary that comes in 3 sizes to help meet the needs of all hospitalized babies. The small internal pellets are smooth and almost feel silky. Unlike more formed boundaries, the Cozy Cub can get into the nooks and crannies to provide support exactly where it is needed. Developed originally in the larger size for ECMO and post-surgical patients, both the large and small Cozy Cub can be used for preemies, agitated babies of all sizes, or any baby who just needs a little 'more.'

SKU	DESCRIPTION	QTY
D 17205 B4	Cozy Cub w/disposable white non-woven cover - Small (40", 3.5 lbs.)	4/case
D 17005 B4	Cozy Cub w/disposable white non-woven cover - Large (48", 6 lbs.)	4/case
D 17805 B2	Cozy Cub w/disposable white non-woven cover - Extra Large (60", 8 lbs.)	2/case
D 17215C B25	Printed cotton covers - Small	25/case
D 17015C B25	Printed cotton covers - Large	25/case
D 17815C B25	Printed cotton covers - Extra Large	25/case
D 17225C B25	Disposable white non-woven covers - Small	25/case
D 17025C B25	Disposable white non-woven covers - Large	25/case
D 17825C B25	Disposable white non-woven covers - Extra Large	25/case

Dandle® PAL, Polly PAL™, & Mini PAL™

Versatile Accessories to Support Neurodevelopmental Positioning



Dandle PAL



Wipeable Dandle



Polly PAL



Mini PAL

Welcome to the PAL family, comprising three flexible positioning aids for NICU babies of all sizes. First is Dandle® PAL, a robust peripheral boundary with a unique design featuring two compartments, or pouches, that allow clinicians to increase height, stability and support where it's needed most. The twill strap stabilizes the smaller pouch on top of or behind the larger pouch, ensuring ample support. Little sister Polly PAL™ features a classic single-pouch design, making it ideal for meeting a myriad of positioning support needs for both the baby and their equipment.

Baby sister Mini PAL™ was welcomed to the family with the unique needs of micro preemies in mind, providing a light weight (only 6oz!) positioning solution that can be used as an insert inside the XS Dandle ROO 2 or ROO Lite, or as a comfort measure with the 'legs' draped gently over the baby. Each member of the PAL family is filled with tiny plastic micro-discs, providing a silky feel. They can be used alone, in pairs, or as an accessory to another positioning aid.

SKU	DESCRIPTION	QTY
D 11505 B10	Dandle PAL- Washable	10/case
D 20505 B10	Dandle PAL- Wipeable, with white disposable cover	10/case
D 20515 C25	Reusable printed cotton covers for Dandle PAL - Wipeable	25/case
D 20525 C25	White disposable covers for Dandle PAL- Wipeable	25/case
D 21755 B10	Polly PAL Infant Positioning Aid	10/case
D 23505 B10	Mini PAL Infant Positioning Aid – Washable	10/case
D 23545 B10	Mini PAL Infant Positioning Aid – Wipeable	10/case
D 23515 B25	Printed Cotton Covers for use with Wipeable Mini PAL	25/case
D 23525 B25	Disposable white non-woven covers for use with Wipeable Mini PAL	25/case

DandleLION® Gel Pillows

A gel pillow that supports developmental positioning

Soft and supportive Gel Pillows in 8 sizes support developmental positioning and help to minimize skin breakdown for the hospitalized baby. All Gel Pillows come with covers. Pillows are not made with natural rubber latex.



SIZE	SKU	DESCRIPTION	QTY
3" x 6"	D 14558 GP B12	Gel Pillow w/ printed cotton covers - 3" x 6"	12/case
	D 14558 GPC B25	Gel Pillow Replacement printed cotton covers - fits 3" x 6" pillow	25/case
	D 12558 GPC B25	Gel Pillow <u>Disposable</u> white covers - fits 3" x 6" pillow	25/case
	D 0958 GP B12	Gel Pillow w/ white cotton covers - 3" x 6" pillow	12/case
6" x 6"	D 0958 GPC B25	Gel Pillow Replacement white cotton covers - fits 3" x 6" pillow	25/case
	D 14559 GP B12	Gel Pillow w/ printed cotton covers - 6" x 6"	12/case
	D 14559 GPC B25	Gel Pillow Replacement printed cotton covers - fits 6" x 6" pillow	25/case
	D 12559 GPC B25	Gel Pillow <u>Disposable</u> white covers - fits 6" x 6" pillow	25/case
6" x 9"	D 0959 GP B12	Gel Pillow w/ white cotton covers - 6" x 6" pillow	12/case
	D 0959 GPC B25	Gel Pillow w/ Replacement white cotton covers - fits 6" x 6" pillow	25/case
	D 14510 GP B12	Gel Pillow w/ printed cotton covers - 6" x 9"	12/case
	D 14510 GPC B25	Gel Pillow Replacement printed cotton covers - fits 6" x 9" pillow	25/case
6" x 11"	D 12551 GPC B25	Gel Pillow <u>Disposable</u> white covers - fits 6" x 9" pillow	25/case
	D 0951 GP B12	Gel Pillow w/ white cotton covers - 6" x 9" pillow	12/case
	D 0951 GPC B25	Gel Pillow Replacement white cotton covers - fits 6" x 9" pillow	25/case
	D 14520 GP B12	Gel Pillow w/ printed cotton covers - 6" x 11"	12/case
6" x 12"	D 14520 GPC B25	Gel Pillow Replacement printed cotton covers - fits 6" x 11" pillow	25/case
	D 12552 GPC B25	Gel Pillow <u>Disposable</u> white covers - fits 6" x 11" pillow	25/case
	D 0952 GP B12	Gel Pillow w/ white cotton covers - 6" x 11"	12/case
	D 0952 GPC B25	Gel Pillow Replacement white cotton covers - fits 6" x 11" pillow	25/case
6" x 7" (~ 3/4" thick)	D 14530 GP B12	Gel Pillow w/ printed cotton cover - 6" x 12"	12/case
	D 14530 GPC B25	Gel Pillow Replacement printed cotton covers - fits 6" x 12" pillow	25/case
	D 12553 GPC B25	Gel Pillow <u>Disposable</u> white covers - fits 6" x 12" pillow	25/case
	D 0953 GP B12	Gel Pillow w/ white cotton cover - 6" x 12"	12/case
6" x 10" (~ 3/4" thick)	D 0953 GPC B25	Gel Pillow Replacement white cotton covers - fits 6" x 12" pillow	25/case
	D 14555 GP B12	Pediatric Gel Pillow w/ printed covers - 6" x 7"	12/case
	D 14555 GPC B25	Pediatric Gel Pillow Replacement printed cotton covers - fits 6" x 7" pillow	25/case
	D 12555 GPC B25	Pediatric Gel Pillow <u>Disposable</u> white covers - fits 6" x 7" pillow	25/case
8" x 13" (~ 1" thick)	D 0955 GP B12	Pediatric Gel Pillow (~7/8" thick) w/ white cotton covers - 6" x 7" pillow	12/case
	D 0955 GPC B25	Pediatric Gel Pillows Replacement white cotton covers - fits 6" x 7" pillow	25/case
	D 14556 GP B12	Pediatric Gel Pillow w/ printed cotton covers - 6" x 10"	12/case
	D 14556 GPC B25	Pediatric Gel Pillow Replacement printed cotton covers - fits 6" x 10" pillow	25/case
8" x 13" (~ 1" thick)	D 12556 GPC B25	Pediatric Gel Pillow <u>Disposable</u> white covers - fits 6" x 10" pillow	25/case
	D 0956 GP B12	Pediatric Gel Pillow w/ white cotton covers - 6" x 10"	12/case
	D 0956 GPC B25	Pediatric Gel Pillow Replacement white cotton covers - fits 6" x 10" pillow	25/case
	D 19054 GP B6	Pediatric Gel Pillow w/ white cotton covers - 8" x 13"	6/case
	D 19054 GPC B25	Pediatric Gel Pillow Replacement white cotton covers - fits 8" x 13" pillow	25/case

Dandy Prone Pad™

Provides Comfortable, Flexed Posture in the Prone Position



The Dandy Prone Pad facilitates developmentally supportive prone positioning that closely mirrors the fetal flexed posture. With its hourglass shape and unique lower extremity supports (goalposts), the Prone Pad encourages the extremities into midline flexion and guides the knees and hips into proper alignment. With 5 sizes that fit easily into DandleLION positioning aids, it's easy to find the right sized support for every patient.

SKU	DESCRIPTION	QTY
D 21531 B12	Prone Pad with 1 Cloth Cover, XS, babies < 1000 g.	12/case
D 21532 B12	Prone Pad with 1 Cloth Cover, S, babies 1000 – 1800 g.	12/case
D 21533 B12	Prone Pad with 1 Cloth Cover, M, babies 1800 – 2500 g.	12/case
D 21534 B12	Prone Pad with 1 Cloth Cover, L, babies 2500 – 3600 g.	12/case
D 21535 B12	Prone Pad with 1 Cloth Cover, XL, babies >3600 g.	12/case
D 21571 B25	Prone Pad Cloth Covers, XS, babies < 1000 g.	25/case
D 21572 B25	Prone Pad Cloth Covers, S, babies 1000 – 1800 g.	25/case
D 21573 B25	Prone Pad Cloth Covers, M, babies 1800 – 2500 g.	25/case
D 21574 B25	Prone Pad Cloth Covers, L, babies 2500 – 3600 g.	25/case
D 21575 B25	Prone Pad Cloth Covers, XL, babies >3600 g.	25/case
D 21591 B50	Prone Pad Disposable Covers, XS, babies < 1000 g.	50/case
D 21592 B50	Prone Pad Disposable Cover, S, babies 1000 – 1800 g.	50/case
D 21593 B502	Prone Pad Disposable Cover, M, babies 1800 – 2500 g.	50/case
D 21594 B50	Prone Pad Disposable Cover, L, babies 2500 – 3600 g.	50/case
D 21595 B50	Prone Pad Disposable Cover, XL, babies >3600 g.	50/case

Organic Skin Care

Cleanse- Nourish- Defend! Organically



Industry-leading formulas made from 70% or more organic ingredients, DandleLION skin care products – Bubbly Wash, Diaper Balm, and Nurturing Oil – are safe and healthy for all babies. They are clinically tested and proven gentle enough for even the most sensitive developing skin. The ingredients are free of dyes, parabens, phthalates, paraffin, silicone, lanolin and irritating essential oils (all the bad stuff!) and each ingredient was chosen for its benefits to the skin.

SKU	DESCRIPTION	QTY
D 10502 Case	BUBBLY WASH - Organic Baby Wash – 10 ml foil packets, Case of 576 (4 boxes of 144)	576/case
D 13645 B40	DIAPER BALM - Organic Diaper Rash Cream - 21 ml flip-top tube, Case of 40 tubes	40/case
D 13620 B96	NURTURING OIL - Organic – 21 ml flip-top tube, Case of 96 tubes	96/case
D 13600 Case	SILKY CREAM - Organic – 10 ml foil packets, Case of 576 (4 boxes of 144)	576/case

DandleLION® Bathing System

Developmental bathing for preemies



A unique bathing system designed specifically for the hospitalized baby . . . the Dandy Tub Complete includes the Dandy Tub with thermometer, a Bathing Hammock that keeps baby flexed and supported during bathing, a disposable liner for easy cleaning and a packet of Organic Bubbly Wash. The Itty Bitty Blanket, that may be purchased as an accessory, helps to swaddle the baby without retaining a lot of water like a washcloth or blanket might do.

SKU	DESCRIPTION	QTY
D 12011 B6	Dandy Tub Complete XSMALL (babies <5 lbs): includes Dandy Tub, 1 XSMALL Bathing Hammock, 1 Plastic Liner, 1-10mL packet of Bubbly Wash	6/case
D 12068 B20	Dandy Tub Wash Packs, XSMALL (babies <5 lbs): includes 1 XSMALL Bathing Hammock, 1 Plastic Liner, 1-10mL packet of Bubbly Wash	20/case
D 12057 B20	Dandy Tub Bathing Hammocks - XSMALL (<5 lbs)	20/case
D 12010 B6	Dandy Tub Complete, SMALL (babies >5 lbs): includes Dandy Tub, 1 SMALL Bathing Hammock, 1 Plastic Liner, 1-10mL packet of Bubbly Wash	6/case
D 12065 B20	Dandy Tub Wash Packs, SMALL (babies >5 lbs): includes 1 SMALL Bathing Hammock, 1 Plastic Liner, 1-10mL packet of Bubbly Wash	20/case
D 12044 B20	Dandy Tub Bathing Hammocks - SMALL (>5 lbs)	20/case
D 12000 B6	Dandy Tub only - no accessories	6/case
D 12033 B100	Dandy Tub plastic liners	100/case
D 12070 B4	Dandy Tub replacement thermometer strips	4/case
D 12155 B25	Itty Bitty Blanket (19" x 21")	25/case

Tranquilo® Soothing Mat

Vibrating Mat to Calm and Soothe Infants



Developed by a labor & delivery nurse, the Tranquilo Soothing Mat mimics mom's sounds and motions to comfort babies who are agitated, stressed or who need sensory input. The vibration levels and pulsatile rhythm choices, which imitate the intrauterine heartbeat, can easily be changed with the push of a button.

SKU	DESCRIPTION	QTY
D 20100 EA	Tranquilo Soothing Mat	1 each
D 20111 B25	Tranquilo Soothing Mat Disposable Covers	25 case

DandleLION® Thermal Hats

A soft warm hug for baby's head



The DandleLION Thermal Hat helps keep newborns from getting cold after delivery, during procedures, transport and in the OR. The quiet Polyethylene Liner (not made with natural rubber latex) is encased between two layers of soft yarn so it never touches the baby. The plastic is stitched securely at the top with overlapping side panels that adjust to the size of the baby's head without compromising the thermal barrier.

SKU	DESCRIPTION	QTY
D 12742 B50	Thermal Hat, Preemie	50/case
D 12843 B50	Thermal Hat, Newborn	50/case
D 17845 B50	Thermal Hat, Newborn Plus	50/case

DandleLION® Sweet Cheeks

40% Glucose Gel




DandleLION is pleased to announce Sweet Cheeks, a 40% glucose gel, made just for babies and designed for clinical convenience. The single-use graduated oral syringe provides controlled administration for all infants, regardless of size. The viscosity of the gel helps it stay on your finger for easy application.

SKU	DESCRIPTION	QTY
D 18200 C200	Sweet Cheeks 40% Glucose Gel - Single Use, - 3mL glucose in 5mL oral syringe- Individually packaged & barcoded (200 syringes)	200/case

DandleLION® KISSES

24% Organic Sucrose Solution



An organic 24% sucrose solution that helps to calm and soothe newborns and babies up to 6 months of age during routine procedures such as blood draws and immunizations. Orange or purple vials: 1 or 2 mL. Eye dropper tip allows controlled administration; each drop is approximately 40-60 microliters.

SKU	DESCRIPTION	QTY
D 12441 C384	DandleLION KISSES (ORANGE) - Single Use, - 1mL vials - Individually bagged & barcoded (384 vials)	384/case
D 12442 C384	DandleLION KISSES (ORANGE) - Single Use - 2mL vials - Individually bagged & barcoded (384 vials)	384/case
D 12481 C384	DandleLION KISSES (PURPLE) - Single Use - 1mL vials - Individually bagged & barcoded (384 vials)	384/case
D 12482 C384	DandleLION KISSES (PURPLE) - Single Use - 2mL vials - Individually bagged & barcoded (384 vials)	384/case



22 Shelter Rock Lane
Danbury, CT 06810
203-791-9001
www.dandlelionmedical.com