

Snuggly ROOTM

Positioning Nest for Increased Visibility

About the Snuggly ROO:

- Quilted base, footroll, and straps made from cotton polyester interlock blend fabric
- · Four sizes:
 - XS: <1000 g.
 - S: 1000-1800 g.
 - · M: 1800-2500 g.
 - · L: 2500-3600 g.
- Each size has a color coded border:
 - XS orange
 - · S yellow
 - M green
 - L blue
- Also available as a disposable made with a soft non-woven material

The Snuggly ROO is part of the DandleLION[®] family of neurodevelopmental positioning products.

It is useful for admissions, transport and the perioperative timeframe, because it provides greater visibility of the baby, as well as support and containment.

The padded footroll keeps the lower extremities flexed and contained, while the offset straps help keep hands to face, shoulders rounded, and arms in flexion. When paired with the Dandle® PAL, the Snuggly ROO is ideal for keeping the head in midline while allowing the baby to remain visible to caregivers. After the early critical period, the baby can transition into the Dandle® ROO2 for more optimal, developmentally appropriate positioning.

Ordering Information

D 21621 B12 | DandleLION Snuggly ROO - XS: <1000 g. | Qty: 12/case
D 21622 B12 | DandleLION Snuggly ROO - S: 1000-1800 g. | Qty: 12/case
D 21623 B12 | DandleLION Snuggly ROO - M: 1800-2500 g. | Qty: 12/case
D 21624 B12 | DandleLION Snuggly ROO - L: 2500-3600 g. | Qty: 12/case
D 21031 B12 | DandleLION Snuggly ROO Disposable - XS: <1000 g. | Qty: 12/case
D 21032 B12 | DandleLION Snuggly ROO Disposable - S: 1000-1800 g. | Qty: 12/case
D 21033 B12 | DandleLION Snuggly ROO Disposable - M: 1800-2500 g. | Qty: 12/case
D 21034 B12 | DandleLION Snuggly ROO Disposable - L: 2500-3600 g. | Qty: 12/case

