



Meets AAP Recommendations for Safe Sleep



Breathable to protect body temperature



Supports babies sleeping on their back



No loose fabric

Dandle® WRAP Stretch

Dynamic Positioning Aid

The Dandle WRAP Stretch features:

- A pouch to maintain lower extremity flexion
- Wide wings that support the shoulders in flexion and encourage hand to mouth self-comforting behaviors
- Velcro® fasteners on pouch and wings that can be adjusted independently of each other to meet baby's individual needs
- Fabric that breathes, yet is absorbent, making it ideal for agitated babies
- Plenty of room between the pouch and the wings for lines or tubing

Babies love to be swaddled! By containing motor activity and mimicking the security of a mother's womb, swaddling comforts babies and results in less crying, more stable heart rates and the ability to sleep through the startle reflex.

In a systematic review of swaddling published in Pediatrics in 2007, authors found "swaddled infants arouse less and sleep longer. Preterm infants have shown improved neuromuscular development, less physiologic distress, better motor organization, and more self-regulatory ability when they are swaddled... (Swaddling) is supportive in cases of neonatal abstinence syndrome and infants with neonatal cerebral lesions."

The Dandle WRAP Stretch offers many swaddling benefits over traditional blankets. Its high performance fabric with 3-way stretch allows freedom of movement while promoting recoil to midline flexion at rest. A wide variety of sizes make this a supportive positioning aid for babies throughout hospitalization.

Swaddling: A Systematic Review van Sleuwen et al, Pediatrics, Vol. 120, No. 4 October 1, 2007, pp. e1097-1106.


Dandle•LION®
MEDICAL

Dandle® WRAP Stretch



Dandle WRAP Stretch is made with fantastic high performance fabric!


- 3 way stretch that returns to its original shape
- Soft, smooth, breathable, wickable, and highly absorbent
- Does not insulate body heat, making it ideal for use in hot and humid environments



Helps during feeding by keeping baby organized and contained. Caregivers can feel respirations through the lightweight fabric and easily see the head and neck in order to visualize sucking and swallowing.



Calms even the most irritable infants. Adjustable swaddling tension and secure Velcro help to calm babies and maintain a flexed and midline position.



Supports during procedures and diaper changes. Wings and pouch can be adjusted independently of each other to keep baby calm and support baby's needs.

D 14531 B12	Dandle WRAP Stretch – XS: <1000 g.	12/case
D 14532 B12	Dandle WRAP Stretch – S: 1000-1800 g.	12/case
D 14533 B12	Dandle WRAP Stretch – M: 1800-2500 g.	12/case
D 14534 B12	Dandle WRAP Stretch – L: 2500-3600 g.	12/case
D 14535 B12	Dandle WRAP Stretch – XL: 3600-4800 g.	12/case
D 14536 B12	Dandle WRAP Stretch – XXL: 4800 g.+	12/case