

Thicker pediatric gel pillows shown with white cotton covers



Standard gel pillows shown with washable printed cotton covers



Dandle•LION® Gel Pillows

With Disposable, White or Printed Cotton Covers

Size	Description
<i>(Standard thickness: ~5/8")</i>	
3" x 6"	Small Roll
6" x 6"	Square Pillow
6" x 9"	Also fits XS Dandle ROO or ROO2
6" x 11"	Also fits S Dandle ROO or ROO2
6" x 12"	Also fits M Dandle ROO or ROO2
<i>(Pediatric or Adult thickness: ~7/8")</i>	
6" x 7"	Pediatric Gel Pillow
6" x 10"	Pediatric Gel Pillow
<i>(Pediatric or Adult thickness: ~1")</i>	
8" x 13"	Pediatric Gel Pillow

Dandle•LION soft gel-filled pillows are designed to provide a comfortable cushion to meet the support needs of micropreemie, preemie, and full term infants.

Gel pillows can help reduce the need for frequent repositioning and make immobilized babies more comfortable. They can also be used as a mattress for very low birthweight babies or as a prone roll for bigger babies.

- Three covers available: white cotton, printed cotton, and disposable non-woven
- Non-toxic, water based food gel
- Guaranteed for six months
- DEHP and BPA free; not made with natural latex rubber
- Single patient use (unless hospital validates their own cleaning process)

Ordering Information

See Product Catalog on resource page of website

