

Polly PAL™

Versatile Positioning Support Filled with Micro Discs

About the Polly PAL:

- · Single bag design
- · Filled with silky micro discs
- Helps achieve and maintain specific positions of baby or equipment
- May be placed around head, neck or hips
- May be partially placed on top of the baby for comfort
- Washable according to CDC guidelines
- Reusable between patients

Caution: Do not place the entire Polly PAL on top of a baby.

Polly PAL, little sister to the Dandle® PAL, is a small, flexible boundary that supports positioning of NICU babies.

The Polly PAL can be used in a myriad of ways, including:

- Around the head to encourage midline positioning
- Supporting respiratory equipment or other medical tubing
- o Behind the back when baby is side-lying
- As a "second set of hands" to provide comfort during care times
- Anywhere the baby needs a little more comfort or support



Ordering Information

D 21755 B10 | Polly PAL Infant Positioning Aid | Qty: 10/case

