

X-SMALL
Less than 6"

SMALL
6" - 7"

MEDIUM
7" - 8.25"

LARGE
8.25" - 9.5"

X-LARGE
More than 9.5"

DANDY PRONE PAD MEASURING GUIDE

Measure the baby from top of head to umbilicus. The approximate length corresponds to the appropriate size.

22 Shelter Rock Lane (p) 203-791-9001
Danbury, CT 06810 (f) 203-792-5581

Dandle·LION
MEDICAL
DandleLIONMedical.com **V2**
Rev3.0424

X-SMALL
Less than 15cm

SMALL
15cm - 17.78cm

MEDIUM
17.78cm - 21cm

LARGE
21cm - 24cm

X-LARGE
More than 24cm

DANDY PRONE PAD MEASURING GUIDE

Measure the baby from top of head to umbilicus. The approximate length corresponds to the appropriate size.

22 Shelter Rock Lane (p) 203-791-9001
Danbury, CT 06810 (f) 203-792-5581

Dandle·LION
MEDICAL
DandleLIONMedical.com **V2**
Rev3.0424