



**Color-coded  
Loop Secures  
Cozy Cub in  
Place**

# Cozy Cub™

## Bead Boundary in 3 Sizes

### Features of the Cozy Cub:

- Easily adaptable to support medical equipment or unique positioning needs
- Capable of providing calming flexion during transport
- Allows immediate access to baby in emergencies
- Sizes:
  - S: 40" long, 3.5 lbs
  - L: 48" long, 6 lbs
  - XL: 60" long, 8 lbs
- Packaged with a disposable cover
- Additional printed cloth covers and disposable white covers available separately

The Cozy Cub is a flexible boundary filled with silky micro discs that nestles comfortably into the nooks and crannies of each baby's unique shape to provide support exactly where it is needed.

It provides complete circumferential support, accommodates a wide variety of medical equipment and allows complete visibility for treatments and procedures such as phototherapy, surgery, or ECMO. The unique securement loop, color-coded by size, maintains the boundary's shape and proximity to the baby, even with movement.

Three sizes allow clinicians to choose the best fit for each baby:

- Small is best suited for ELBW, VLBW and SGA babies
- Large is best suited for preterm and term newborns of all gestational ages
- Extra large is best suited for term newborns and older babies

### Ordering Information

D 17205 B4 | Cozy Cub and white non-woven cover - Small | Qty: 4/case

D 17005 B4 | Cozy Cub and white non-woven cover - Large | Qty: 4/case

D 17805 B4 | Cozy Cub and white non-woven cover - XL | Qty: 2/case

See website for part numbers for additional covers, including washable cloth or non-woven disposable styles.

