

Patented PneuMaPure® filter technology allows airflow while blocking pathogens such as MRSA, C. diff, E. coli, coronavirus, and more



Pureshield™ fabric creates a waterproof exterior preventing pathogen ingress and makes mattress easy-to-clean with no laundry needed



Safeweld® hermetically-sealed seams keep pathogens out by eliminating stitching holes

Dandle® PURE Neonatal Mattress

With PneuMaPure® filter

DandlePURE Neonatal Mattress features:

- Fits standard radiant warmers, incubators, bassinets, and cribs
- Multi-patient use
- Pureshield fabric is biocompatible
- Fabric and integrated filter can be wiped clean with approved hospital disinfectants
- Core made of viscoelastic foam (memory foam)
- DEHP and latex-free
- 2-year warranty against manufacturer defects

A cleaner, safer alternative to traditional NICU mattresses.

Rest easy, knowing your babies are safer from the start. DandlePURE Neonatal Mattresses with PneuMaPure filter technology can replace your standard mattress while offering 3-way protection from harmful pathogens. PneuMaPure is a clinically-proven, lab-tested filter for soft surfaces.


DandlePURE protects patients through:

- An integrated PneuMaPure filter that prevents ingress of pathogens and allergens while allowing air to escape, keeping the mattress soft and comfortable.
- Soft, durable, high-performance Pureshield fabric not only blocks out contaminants, but is waterproof.
- Safeweld air-tight seams (no stitching holes), which are heat-sealed to provide a complete barrier.

Rest assured. DandlePURE.



Dandle®PURE Neonatal Mattress



The patented PneuMapure® filter technology blocks bacteria, allergens, mold, and deadly viruses while maintaining sufficient breathability.

| | | |
|---------|--|------|
| D 24770 | Dandle®PURE Neonatal Mattress with PneuMapure® Filter: fits Panda™ Warmer | 1 EA |
| D 24771 | Dandle®PURE Neonatal Mattress with PneuMapure® Filter: fits Panda™ L&D Warmer | 1 EA |
| D 24780 | Dandle®PURE Neonatal Mattress with PneuMapure® Filter: fits Giraffe™ Family | 1 EA |
| D 24790 | Dandle®PURE Neonatal Mattress with PneuMapure® Filter: fits Airborne® 750i | 1 EA |
| D 24800 | Dandle®PURE Neonatal Mattress with PneuMapure® Filter: fits Bassinet | 1 EA |
| D 24810 | Dandle®PURE Neonatal Mattress with PneuMapure® Filter: fits Cozy Cot Warmer | 1 EA |
| D 24820 | Dandle®PURE Neonatal Mattress with PneuMapure® Filter: fits Air Shields® Isolette® | 1 EA |
| D 24830 | Dandle®PURE Neonatal Mattress with PneuMapure® Filter: fits ATOM Dual Incu i | 1 EA |
| D 24840 | Dandle®PURE Neonatal Mattress with PneuMapure® Filter: fits NÄRA Bassinet | 1 EA |
| D 24850 | Dandle®PURE Neonatal Mattress with PneuMapure® Filter: fits Caleo® Dräger | 1 EA |
| D 24860 | Dandle®PURE Neonatal Mattress with PneuMapure® Filter: fits Hard Crib Large | 1 EA |
| D 24870 | Dandle®PURE Neonatal Mattress with PneuMapure® Filter: fits Hard Crib Small | 1 EA |

