

- Ammonium Compounds
- Aniospray Quick
- AntiGone
- Bacillol AF
- Bacillol Foam
- Bactacid Wipes
- Chemisept Wipes
- Chlorine Dioxide Solution, 10,000 ppm (for use in the case of heavy soiling; clean with a damp cloth after use)
- Chlorine Dioxide Solution, 1000 ppm
- Chlorine Releasing Solution, 1000 ppm
- Clinell Alcohol Wipes
- Clinell Sporidical Wipe
- Clinell Universal Spray
- Cavi Cide Surface Decontaminant Cleaner
- Cavi Wipes
- Clorox Bleach Germicidal Wipes, NaOCl 0.55%
- Descosept PUR Wipes RTU
- Ecopine, (5% Solution in Water)
- Hydrogen Peroxide Solution, 11%
- Incidin Foam & Incidin Plus
- Incidin OxyFoam S
- Incidin OxyWipes
- Incidin Pro
- Isopropyl Alcohol, 70%
- Kohrsolin FF
- Life Clean
- Meliseptol Rapid
- Microbac Forte
- Mikrozyd AF Liquid
- Mikrozyd Sensitive Liquid
- Mikrozyd Sensitive Wipes
- Oxivir 1
- Oxivir Plus
- Oxivir TB
- Oxycide
- Perform 1% and 3%
- Phenol Solution, 1%
- PDI Sani Cloth AF3 Wipes
- PDI Sani-Cloth Bleach Germicidal Disposable Wipe
- Rely+ON Perasafe
- Rely+ON Virkon
- Sani-Cloth Active
- Spray Off
- Sterisept Wipes
- Super Sani Cloth Germicidal Disposable Wipe
- Terralin® Protect (max. 4%)
- TPH Protect
- Virex 256
- Virusolve (5% & 10% Solution in Water) Wet Wipe

Note: Chlorine can damage the surface fabric of the product during long-term use and make it more sensitive to mechanical damage. After using disinfectants containing chlorine, it is recommended to wipe the surfaces with a damp cloth.

Not recommended: Disclean ENDO PAA