



Hands Free
Collection Cups
available!



Ardo Bellis Breast Pump

Hospital-Grade Multi-User Personal Breast Pump

Key Features:

- Ideal for use in a single hospital room or as part of a rental fleet
- Whisper-quiet motor (~40 db) provides
 - A max suction of 250 mmHg
 - A max speed of 110 cycles per minute
- Completely closed system
- Automatic Power Pumping mode
- Memory Plus functionality
- Seamless exterior for effective cleaning between patients
- Rechargeable integrated battery allows portability
- Compact (10"x8"x6")
- Weighs 1.4 pounds
- Built-in bottle holder/handle
- Works as a single or double pump
- Compatible with Ardo Hands-Free Collection Cups
- Universal pump set connects to any Ardo pump
- Wide variety of flange sizing options
- BPA & BPS free
- FDA-cleared
- Connects to MyArdo smartphone app

Meet Ardo Bellis — the world's first and only multi-user breast pump to feature an Automatic Power Pumping function that mimics cluster feeding, which is shown to increase milk supply. This hospital-grade, multi-user breast pump is tailored for mothers facing the challenges of nurturing hospitalized infants.

Weighing just 1.4 pounds and featuring a rechargeable battery, its compact and portable design ensures pumping convenience anywhere without compromising pumping power or milk supply. This makes Bellis the ideal take-home companion for your rental fleet.


Mothers can operate Bellis using the free MyArdo app and save their preferences with the Memory Plus function in the app or on the pump face.

The completely closed system, with the barrier built into the pump set, safeguards against milk entering the tubing or pump. Coupled with a sleek touchscreen panel for efficient cleaning, the Bellis sets a new hygienic standard, surpassing the expectations of today's moms and hospital teams alike.

Dandle•LION[®]
MEDICAL

Distributed by Dandle•LION Medical

Ardo Bellis Breast Pump



MyArdo App Features:

- Start/stop the pump via smartphone
- All information about operating the breast pump
- Operation of the Power Pumping function
- Operation of the Memory Plus function
- Vacuum and cycle setting
- Tips and tricks on pumping

ARDO 63.00.524	Ardo Bellis hospital-grade multi-user Swiss made breast Pump. For hospital & home use by one mom at a time. Ideal for private rooms, rental or loaner programs.	1/box
ARDO 63.00.114 Q30	One-Mom Pump Set D26 - 26mm shell, milk bottle 150ml (5 fl.oz): Complete with valve, tube, membrane & cap, Sterile (1 set, 1 pack)	30 sets/box
ARDO 63.00.125 Q30	One-Mom Pump Set D31 - 31mm shell, milk bottle 150ml (5 fl.oz): Complete with valve, tube, membrane & cap, Sterile (1 set, 1 pack)	30 sets/box
ARDO 63.00.239 Q25	Amaryll Essentials Hand Pump - 26mm shell, 22mm hard insert, milk bottle 150mL, Complete with lever, valve, membrane & rod, cap & stand, non-sterile (25 units)	25/box
ARDO 63.00.454 Q10	Hands-Free Pump Set - 2 wearable milk collection cup sets, fits all Ardo electric breast pumps (non-sterile)	10 sets/box
ARDO 63.00.555 Q45	Double Clinic Pumpset; 26mm shell, milk bottles 150ml (5 fl.oz): Complete with valves, tubes, membranes & caps, Non-sterile (2 sets, 1 pack)	45 sets/box