



Hands-Free  
Collection Cups  
available!



# Ardo Carum Breast Pump

## Hospital-Grade Breast Pump

### Key Features:

- Whisper-quiet dual piston motor provides
  - A max suction of 250 mmHg (while double pumping)
  - A max speed of 120 cycles per minute
- Completely closed system
- Sensitive Programme™
- Stimulation and expression modes for 2-phase function
- Large, color LCD display
- Seamless exterior for effective cleaning between patients
- Built-in bottle holder
- Mounts to sturdy rolling stand or sits on a countertop
- Works as a single or double pump
- Compatible with Ardo Hands-Free Collection Cups
- Universal pump set connects to any Ardo pump
- Wide variety of flange sizing options
- BPA & BPS free
- FDA-cleared

Meet Ardo Carum, an innovative hospital-grade breast pump that caters to the needs of today's moms while meeting the standards of modern hospitals.

Mothers of hospitalized infants can pump confidently and comfortably with the Carum: a quiet, customizable, and hygienic upgrade to standard multi-user, hospital-grade breast pumps.

Its powerful dual-piston motor is remarkably quiet, minimizing ambient noise, making it the ideal pump for quiet hospital settings. Customizable features, including the Sensitive Programme™ and the Let-Down Button, cater to individual needs, particularly for new moms or those with premature births. With Carum, mothers can select suction and speed independently from each other, unlike other leading hospital pumps.

Ardo Carum's completely closed system with the barrier built into the pump set ensures milk stays out of the tubing and pump. Complemented by its sleek exterior, Carum presents an intuitively safer and more hygienic solution.

Revolutionize hospital-grade pumping with the Ardo Carum—for superior infection control, customization, and quiet efficiency, benefitting both clinicians and mothers alike.

**Dandle•LION<sup>®</sup>**  
MEDICAL

Distributed by Dandle•LION Medical

# Ardo Carum Breast Pump



ARDO 63.00.73	Ardo Carum hospital-grade multi-user breast pump made in Switzerland. For 24/7 hospital operation w/ multiple mothers.	1/box
ARDO 63.00.129	Mobile Stand for Ardo Carum hospital pump: includes two bottle holders	1/box
ARDO 63.00.114 Q30	One-Mom Pump Set D26 - 26mm shell, milk bottle 150ml (5 fl.oz): Complete with valve, tube, membrane & cap, Sterile (1 set, 1 pack)	30 sets/box
ARDO 63.00.125 Q30	One-Mom Pump Set D31 - 31mm shell, milk bottle 150ml (5 fl.oz): Complete with valve, tube, membrane & cap, Sterile (1 set, 1 pack)	30 sets/box
ARDO 63.00.239 Q25	Amaryll Essentials Hand Pump - 26mm shell, 22mm hard insert, milk bottle 150mL, Complete with lever, valve, membrane & rod, cap & stand, non-sterile (25 units)	25/box
ARDO 63.00.454 Q10	Hands-Free Pump Set - 2 wearable milk collection cup sets, fits all Ardo electric breast pumps (non-sterile)	10 sets/box
ARDO 63.00.555 Q45	Double Clinic Pumpset; 26mm shell, milk bottles 150ml (5 fl.oz): Complete with valves, tubes, membranes & caps, Non-sterile (2 sets, 1 pack)	45 sets/box